

A HYGIENIC PARADIGM SHIFT



WISHOWER

By LIOMAR Medical



AGENDA

- Other product lines
- About LIOMAR Medical
- Challenges and risks of traditional hygienic care for patients in care units
- Wishower
 - A hygienic paradigm shift
 - Humanization of hygienic care
 - Types of patients
 - Components and accessories
 - Wishower video
 - Features and benefits
 - Disinfection process
- Questions/Answers

We also represent...



Wishower Shower in bed

Mobile shower in bed solution
mobile for mobility reduced and
bedridden patients.



SafeCrush automatic pill crusher

SafeCrush stands out with its
automatic, efficient, safe and
silent crushing system.



Wide range of PRIMED Medical PPE

Many innovative and industry-
leading Canadian product lines
used in many Quebec facilities.



ABOUT US

Products and positioning

New innovative technologies based on a deep understanding of both patients in care units and healthcare workers.

Company background

Canadian based company led by executives cumulating many years of experience in the healthcare industry.

Mission

Improving vulnerable patients lives

Market

Wishower is currently being used in many healthcare facilities across Canada (Quebec, Ontario and Atlantic region). It is now also available in the U.S.

LIOMAR
M E D I C A L



CHALLENGES AND RISKS

Associated with traditional hygienic care for patients in care units and for healthcare workers

CHALLENGES

Patients requiring special assistance

Concerns related to the medical conditions of patients

Obsolete or inadequate equipment

Complex or difficult equipment to operate

Lack of nursing staff

Direct impact on the quality care provided

Pandemic

Active outbreaks, limited resources and time allowed to the disinfection

Requirements and expectations

Of patients, their families, government standards, etc.

RISKS



Handling and maneuvering

Discomfort, risk of accidents, risk of injuries and compromised efficiency

Infection spread

Inter-patient transport complex protocols and cross-contamination

When patients are moved to the bathing area, there's an inherent risk of both infection transmission and injuries.

Inadequate for some patients

The experience can be traumatic and impact patient's dignity

THE WISHOWER SOLUTION

A paradigm shift in hygienic care for
vulnerable patients and healthcare
workers



Health Canada
Certified



Compliant with
CSA standards



Designed and
made in Canada

HYGIENIC PARADIGM SHIFT



« We do not need to move the patient to the hygienic care bathing area. We bring the hygienic care solution to the patient. »

HUMANIZATION OF HYGIENIC CARE

OUR **THREE** PILLARS



Dignity

The patient experiences a complete shower in the intimacy and comfort of their own bed.



Safety


Mobile solution that minimizes the risks of accidents, injuries and contamination.



Ergonomy

The Wishower requires minimum effort to operate and move around the facility.

TYPES OF PATIENTS




Who can benefit from Wishower?

This device can provide complete hygienic care to all patients needing special assistance.

A selection of patients that are more likely to benefit from this type of care is recommended.

DEVICE COMPONENTS



WISHOWER VIDEO



WISHOWER ACCESSORIES



Reusable
Bed size towel



Reusable or **disposable**
waterproof bed cover



Reusable Set of
four (4) clamps



OPTIONAL ACCESSORIES



Disposable dry wash glove



Reusable shampoo basin

FEATURES AND BENEFITS

Patient's well-being and comfort

- Adjustable water temperature according to patient preference (37°C-38°C)
- Constant water temperature control
- Consistant water jet pressure
- Low noise sound profile (30 dBA)
- Performant suction tip designed to only suction liquids and residues

Simple and easy use for staff

- User friendly tactile interface that displays instructions and visual indicators (EN, FR, Spanish)
- Simultaneous water spray/suction
- Integrated and automated disinfection system
- One full water tank of **28L for up to 2 patients**
- Complete care including disinfection process can be done by one person in **under 30 minutes**
- Complete in-bed hygienic care takes **8 minutes** per patient

The Wishower comes with a **2-year warranty**

DISINFECTION PROCESS

Integrated system

1. Self contained and automated cycle
2. Process takes 4-5 minutes
3. Step by step instructions through the disinfection process

Outer surfaces

1. Smooth and easy to clean surfaces with disinfectant wipes
2. Components and accessories easy to clean and disinfect (hoses, clamps, removable filter)
3. Can be done during the internal cycle

Disinfecting solutions approved by healthcare facilities can be used for both the **inner integrated system** and **outer surfaces** (ex : Oxyvir, hydrogen peroxyde).

DISINFECTION VIDEO



PRODUCTIVITY AND EFFICIENCY COMPARISON

	WISHOWER	STRETCHER SHOWER/SHOWE R GURNEY	SHOWER	THERAPEAUTIC BATH
Water consumption	14L	60L	75L	160L
Complete process (shower and disinfection)	+/- 26 min	+/- 45 min	+/- 40 min	+/- 45 min
Personnel	1-2 staff member	* 2-3 staff members	* 2-3 staff members	* 2-3 staff members

* Based on 2 healthcare workers prepping/tranportation/bathing and
1 maintenance person for the disinfection process.

Source : Compiled data from internal tests with multiple Canadian LTC/nursing homes.

Quick reminder...



Wishower Shower in bed

Mobile shower in bed solution
mobile for mobility reduced and
bedridden patients.



SafeCrush automatic pill crusher

SafeCrush stands out with its
automatic, efficient, safe and
silent crushing system.



Wide range of PRIMED Medical PPE

Many innovative and industry-
leading Canadian product lines
used in many Quebec facilities.



THE HUMANIZATION OF HYGIENIC CARE

LIOMAR
MEDICAL

Your Wishower Representatives from LIOMAR Medical

Diane Sawaya

Vice-President, Associate

✉ diane@liomarmedical.com

☎ 514-770-7323

Camille Bergeron

Assistant Manager Business development

✉ camille@liomarmedical.com

Sherry Chau-Carganilla

Customer Service Representative

✉ sherry.chau@liomarmedical.com



www.wishower.net



+1 450-987-1127



info@liomarmedical.co